

Mediator of RNA polymerase II transcription subunit 31

Catalog No: #AP85669

Orders: order@signalwayantibody.com

Package Size: #AP85669-1 50ug #AP85669-2 100ug #AP85669-3 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Mediator of RNA polymerase II transcription subunit 31
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Mediator complex subunit 31, Mediator complex subunit SOH1, hSOH1, SOH1, CGI-125
Accession No.	Uniprot ID:Q9Y3C7
Uniprot	Q9Y3C7
GeneID	51003;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: MED31

Note: This product is for in vitro research use only