

S-adenosylmethionine synthase isoform type-1

Catalog No: #AP85670



Package Size: #AP85670-1 50ug #AP85670-2 100ug #AP85670-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	S-adenosylmethionine synthase isoform type-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	AdoMet synthase 1,Methionine adenosyltransferase 1,MAT 1,Methionine adenosyltransferase I/III,MAT-I/III,Ams1
Accession No.	Uniprot ID:P13444
Uniprot	P13444
GeneID	25331;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mat1a

Note: This product is for in vitro research use only