

N-alpha-acetyltransferase 20

Catalog No: #AP85672



Package Size: #AP85672-1 50ug #AP85672-2 100ug #AP85672-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	N-alpha-acetyltransferase 20
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Methionine N-acetyltransferase,N-acetyltransferase 5,N-terminal acetyltransferase B complex catalytic subunit NAA20,N-terminal acetyltransferase B complex catalytic subunit NAT5,NatB complex subunit NAT5,NatB catalytic subunit,Nat5
Accession No.	Uniprot ID:P61600
Uniprot	P61600
GeneID	67877;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Naa20

Note: This product is for in vitro research use only