

NMDA receptor-regulated protein 2

Catalog No: #AP85674



Package Size: #AP85674-1 50ug #AP85674-2 100ug #AP85674-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	NMDA receptor-regulated protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Accession No.	Uniprot ID:Q0VCQ7
Uniprot	Q0VCQ7
GeneID	541255;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NARG2

Note: This product is for in vitro research use only