

Dual specificity mitogen-activated protein kinase kinase 1

Catalog No: #AP85675

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP85675-1 50ug #AP85675-2 100ug #AP85675-3 1mg

Description

Product Name	Dual specificity mitogen-activated protein kinase kinase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	MAP kinase kinase 1,ERK activator kinase 1,MAPK/ERK kinase 1,MEK 1,Mek1,Prkmk1
Accession No.	Uniprot ID:Q01986
Uniprot	Q01986
GeneID	170851;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Map2k1

Note: This product is for in vitro research use only