

## Mucin-1

Catalog No: #AP85676



Package Size: #AP85676-1 50ug #AP85676-2 100ug #AP85676-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mucin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MUC-1,Breast carcinoma-associated antigen DF3,Cancer antigen 15-3,CA 15-3,Carcinoma-associated mucin,Episialin,H23AG,Krebs von den Lungen-6,KL-6,PEMT,Peanut-reactive urinary mucin,PUM,Polymorphic epithelial mucin,PEM,Tumor-associated epithelial membrane antigen,EMA,Tumor-associated mucin,CD227,PUM
Accession No.	Uniprot ID:P15941
Uniprot	P15941
GeneID	4582;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:MUC1

Note: This product is for in vitro research use only