

Mucin-15

Catalog No: #AP85677



Package Size: #AP85677-1 50ug #AP85677-2 100ug #AP85677-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mucin-15
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	MUC-15,Component II,Glycoprotein 4,Glycoprotein C,PAS3,PASIII
Accession No.	Uniprot ID:Q8MI01
Uniprot	Q8MI01
GeneID	337919;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MUC15

Note: This product is for in vitro research use only