

Matrix metalloproteinase-9

Catalog No: #AP85680



Package Size: #AP85680-1 50ug #AP85680-2 100ug #AP85680-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Matrix metalloproteinase-9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	MMP-9,92 kDa gelatinase,92 kDa type IV collagenase,Gelatinase B,GELB
Accession No.	Uniprot ID:P52176
Uniprot	P52176
GeneID	282871;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MMP9

Note: This product is for in vitro research use only