

Nuclear receptor coactivator 1

Catalog No: #AP85682



Package Size: #AP85682-1 50ug #AP85682-2 100ug #AP85682-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nuclear receptor coactivator 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	NCoA-1,Nuclear receptor coactivator protein 1,mNRC-1,Steroid receptor coactivator 1,Src1
Accession No.	Uniprot ID:P70365
Uniprot	P70365
GeneID	17977;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ncoa1

Note: This product is for in vitro research use only