

Lutropin subunit beta

Catalog No: #AP85684



Package Size: #AP85684-1 50ug #AP85684-2 100ug #AP85684-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Lutropin subunit beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Luteinizing hormone subunit beta,LH-B,LSH-B,LSH-beta,Lutropin beta chain
Accession No.	Uniprot ID:P18842
Uniprot	P18842
GeneID	403959;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LHB

Note: This product is for in vitro research use only