

Major vault protein

Catalog No: #AP85685



Package Size: #AP85685-1 50ug #AP85685-2 100ug #AP85685-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Major vault protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MVP,Lung resistance-related protein,LRP
Accession No.	Uniprot ID:Q14764
Uniprot	Q14764
GeneID	9961;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MVP

Note: This product is for in vitro research use only