

Myogenic factor 6

Catalog No: #AP85686



Package Size: #AP85686-1 50ug #AP85686-2 100ug #AP85686-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Myogenic factor 6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Myf-6,Herculin,Muscle-specific regulatory factor 4,Mrf4,Myf-6
Accession No.	Uniprot ID:P15375
Uniprot	P15375
GeneID	17878;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Myf6

Note: This product is for in vitro research use only