

NAD(P)H-hydrate epimerase

Catalog No: #AP85693



Package Size: #AP85693-1 50ug #AP85693-2 100ug #AP85693-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | NAD(P)H-hydrate epimerase |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Dog |
| Other Names | Apolipoprotein A-I-binding protein,AI-BP,NAD(P)HX epimerase,AIBP |
| Accession No. | Uniprot ID:E2QRY6 |
| Uniprot | E2QRY6 |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:APOA1BP

Note: This product is for in vitro research use only