

Mesoderm-specific transcript homolog protein

Catalog No: #AP85698



Package Size: #AP85698-1 50ug #AP85698-2 100ug #AP85698-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mesoderm-specific transcript homolog protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Paternally-expressed gene 1 protein, PEG1
Accession No.	Uniprot ID:Q2HJM9
Uniprot	Q2HJM9
GeneID	404180;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: MEST

Note: This product is for in vitro research use only