

Neurogenic locus notch homolog protein 2

Catalog No: #AP85706



Package Size: #AP85706-1 50ug #AP85706-2 100ug #AP85706-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neurogenic locus notch homolog protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Notch 2
Accession No.	Uniprot ID:Q9QW30
Uniprot	Q9QW30
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Notch2

Note: This product is for in vitro research use only