

# Mucin-2

Catalog No: #AP85707



Package Size: #AP85707-1 50ug #AP85707-2 100ug #AP85707-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Mucin-2  |
| Brief Description  | Recombinant Protein                                    |
| Host Species       | E.coli   |
| Purification       | Greater than 90% by SDS-PAGE                           |
| Species Reactivity | Ms   |
| Other Names        | MUC-2,Colonic mucin,MCM,Secreted gel-forming mucin     |
| Accession No.      | Uniprot ID:Q80Z19                                      |
| Uniprot            | Q80Z19   |
| Tag Info           | His  |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.) |

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Muc2

Note: This product is for in vitro research use only