

# Microtubule-associated protein 2

Catalog No: #AP85708



Package Size: #AP85708-1 50ug #AP85708-2 100ug #AP85708-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Microtubule-associated protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	MAP-2,Mtap2
Accession No.	Uniprot ID:P20357
Uniprot	P20357
GeneID	17756;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Map2

Note: This product is for in vitro research use only