

Neural cell adhesion molecule 1

Catalog No: #AP85711



Package Size: #AP85711-1 50ug #AP85711-2 100ug #AP85711-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neural cell adhesion molecule 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	N-CAM-1,NCAM
Accession No.	Uniprot ID:P31836
Uniprot	P31836
GeneID	281941;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NCAM1

Note: This product is for in vitro research use only