

Protein NDRG1

Catalog No: #AP85712



Package Size: #AP85712-1 50ug #AP85712-2 100ug #AP85712-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein NDRG1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	N-myc downstream-regulated gene 1 protein, Protein Ndr1, Ndr1, Ndr1, Tdd5
Accession No.	Uniprot ID:Q62433
Uniprot	Q62433
GeneID	17988;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Ndr1

Note: This product is for in vitro research use only