

2'-5'-oligoadenylate synthase 1

Catalog No: #AP85714



Package Size: #AP85714-1 50ug #AP85714-2 100ug #AP85714-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	2'-5'-oligoadenylate synthase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	(2-5')oligo(A) synthase 1,p42 OAS
Accession No.	Uniprot ID:Q29599
Uniprot	Q29599
GeneID	397570;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: OAS1

Note: This product is for in vitro research use only