

## Phosphate carrier protein, mitochondrial

Catalog No: #AP85716



Package Size: #AP85716-1 50ug #AP85716-2 100ug #AP85716-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Phosphate carrier protein, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Phosphate transport protein,PTP,Solute carrier family 25 member 3
Accession No.	Uniprot ID:Q8VEM8
Uniprot	Q8VEM8
GeneID	18674;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Slc25a3

Note: This product is for in vitro research use only