

Neudesin

Catalog No: #AP85720



Package Size: #AP85720-1 50ug #AP85720-2 100ug #AP85720-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neudesin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Neuron-derived neurotrophic factor, SCIRP10-related protein, Spinal cord injury-related protein 10
Accession No.	Uniprot ID:Q1JQA5
Uniprot	Q1JQA5
GeneID	616334;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: NENF

Note: This product is for in vitro research use only