

# Mitogen-activated protein kinase 1

Catalog No: #AP85723



Package Size: #AP85723-1 50ug #AP85723-2 100ug #AP85723-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Mitogen-activated protein kinase 1   |
| Brief Description  | Recombinant Protein  |
| Host Species       | E.coli   |
| Purification       | Greater than 90% by SDS-PAGE   |
| Species Reactivity | Hu   |
| Other Names        | MAP kinase 1,ERT1,Extracellular signal-regulated kinase 2,ERK-2,MAP kinase isoform p42,p42-MAPK,Mitogen-activated protein kinase 2,MAP kinase 2,ERK2,PRKM1,PRKM2 |
| Accession No.      | Uniprot ID:P28482  |
| Uniprot            | P28482   |
| GeneID             | 5594;  |
| Tag Info           | His  |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)   |

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:MAPK1

Note: This product is for in vitro research use only