

Mitogen-activated protein kinase 15

Catalog No: #AP85724



Package Size: #AP85724-1 50ug #AP85724-2 100ug #AP85724-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mitogen-activated protein kinase 15
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MAP kinase 15, Extracellular signal-regulated kinase 7, ERK-7, Extracellular signal-regulated kinase 8, ERK-8, ERK7, ERK8
Accession No.	Uniprot ID: Q8TD08
Uniprot	Q8TD08
GeneID	225689;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: MAPK15

Note: This product is for in vitro research use only