

Matrix metalloproteinase-24

Catalog No: #AP85726



Package Size: #AP85726-1 50ug #AP85726-2 100ug #AP85726-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Matrix metalloproteinase-24
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MMP-24,Membrane-type matrix metalloproteinase 5,MT-MMP 5,Membrane-type-5 matrix metalloproteinase,MT5-MMP,MT5MMP
Accession No.	Uniprot ID:Q9Y5R2
Uniprot	Q9Y5R2
GeneID	10893;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MMP24

Note: This product is for in vitro research use only