

# Dual specificity mitogen-activated protein kinase kinase 6

Catalog No: #AP85727

Package Size: #AP85727-1 50ug #AP85727-2 100ug #AP85727-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Dual specificity mitogen-activated protein kinase kinase 6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MAP kinase kinase 6,MAPK/ERK kinase 6,MEK 6,Stress-activated protein kinase kinase 3,SAPK kinase 3,MEK6,MKK6,PRKMK6,SKK3
Accession No.	Uniprot ID:P52564
Uniprot	P52564
GeneID	5608;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:MAP2K6

Note: This product is for in vitro research use only