

Neuroigin-3

Catalog No: #AP85728



Package Size: #AP85728-1 50ug #AP85728-2 100ug #AP85728-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Neuroigin-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Gliotactin homolog,KIAA1480,NL3
Accession No.	Uniprot ID:Q9NZ94
Uniprot	Q9NZ94
GeneID	54413;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NLGN3

Note: This product is for in vitro research use only