

## Neurogenic locus notch homolog protein 1

Catalog No: #AP85729



Package Size: #AP85729-1 50ug #AP85729-2 100ug #AP85729-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Neurogenic locus notch homolog protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Notch 1, Translocation-associated notch protein TAN-1, TAN1
Accession No.	Uniprot ID:P46531
Uniprot	P46531
GeneID	4851;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: NOTCH1

Note: This product is for in vitro research use only