

Mucin-20

Catalog No: #AP85731



Package Size: #AP85731-1 50ug #AP85731-2 100ug #AP85731-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mucin-20
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MUC-20,KIAA1359,UNQ2782/PRO7170
Accession No.	Uniprot ID:Q8N307
Uniprot	Q8N307
GeneID	200958;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MUC20

Note: This product is for in vitro research use only