

NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 6

Catalog No: #AP85732

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP85732-1 50ug #AP85732-2 100ug #AP85732-3 1mg

Description

Product Name	NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Complex I-B14,CI-B14,LYR motif-containing protein 6,NADH-ubiquinone oxidoreductase B14 subunit,LYRM6,NADHB14
Accession No.	Uniprot ID:P56556
Uniprot	P56556
GeneID	4700;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NDUFA6

Note: This product is for in vitro research use only