

Myelin-associated glycoprotein

Catalog No: #AP85734



Package Size: #AP85734-1 50ug #AP85734-2 100ug #AP85734-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Myelin-associated glycoprotein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	1B236,Brain neuron cytoplasmic protein 3,Sialic acid-binding Ig-like lectin 4a,Siglec-4a
Accession No.	Uniprot ID:P07722
Uniprot	P07722
GeneID	29409;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mag

Note: This product is for in vitro research use only