

# Myosin-11

Catalog No: #AP85735



Package Size: #AP85735-1 50ug #AP85735-2 100ug #AP85735-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Myosin-11
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Myosin heavy chain 11, Myosin heavy chain, smooth muscle isoform, SMMHC
Accession No.	Uniprot ID:P81271
Uniprot	P81271
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: MYH11

Note: This product is for in vitro research use only