

Melanopsin

Catalog No: #AP85736



Package Size: #AP85736-1 50ug #AP85736-2 100ug #AP85736-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Melanopsin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Opsin-4,Mop,Mopn
Accession No.	Uniprot ID:Q9QXZ9
Uniprot	Q9QXZ9
GeneID	30044;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Opn4

Note: This product is for in vitro research use only