

Dehydrodolichyl diphosphate synthase complex subunit NUS1

Catalog No: #AP85738

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP85738-1 50ug #AP85738-2 100ug #AP85738-3 1mg

Description

Product Name	Dehydrodolichyl diphosphate synthase complex subunit NUS1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Di-trans,poly-cis-decaprenylcistransferase,Nogo-B receptor,NgBR,Nuclear undecaprenyl pyrophosphate synthase 1 homolog,C6orf68,NGBR
Accession No.	Uniprot ID:Q96E22
Uniprot	Q96E22
GeneID	116150;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NUS1

Note: This product is for in vitro research use only