

Hepatocyte growth factor receptor

Catalog No: #AP85747



Package Size: #AP85747-1 50ug #AP85747-2 100ug #AP85747-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Hepatocyte growth factor receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	HGF receptor,HGF/SF receptor,Proto-oncogene c-Met,Scatter factor receptor,SF receptor,Tyrosine-protein kinase Met
Accession No.	Uniprot ID:Q75ZY9
Uniprot	Q75ZY9
GeneID	403438;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MET

Note: This product is for in vitro research use only