

Lipopolysaccharide-induced tumor necrosis factor-alpha factor

Catalog No: #AP85748

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP85748-1 50ug #AP85748-2 100ug #AP85748-3 1mg

Description

Product Name	Lipopolysaccharide-induced tumor necrosis factor-alpha factor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	LPS-induced TNF-alpha factor, Small integral membrane protein of lysosome/late endosome, p53-induced gene 7 protein, PIG7, SIMPLE
Accession No.	Uniprot ID: Q99732
Uniprot	Q99732
GeneID	9516;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: LITAF

Note: This product is for in vitro research use only