

# Metallo-beta-lactamase domain-containing protein 1

Catalog No: #AP85749



Package Size: #AP85749-1 50ug #AP85749-2 100ug #AP85749-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Metallo-beta-lactamase domain-containing protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Accession No.	Uniprot ID:A4D2B0
Uniprot	A4D2B0
GeneID	255374;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:MBLAC1

Note: This product is for in vitro research use only