

Membrane-bound transcription factor site-1 protease

Catalog No: #AP85750



Package Size: #AP85750-1 50ug #AP85750-2 100ug #AP85750-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Membrane-bound transcription factor site-1 protease
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Endopeptidase S1P, Sterol-regulated luminal protease, Subtilisin/kexin isozyme 1, SKI-1, S1p, Ski1
Accession No.	Uniprot ID:Q9WTZ2
Uniprot	Q9WTZ2
GeneID	56453;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Mbtps1

Note: This product is for in vitro research use only