

Melanocortin receptor 5

Catalog No: #AP85751



Package Size: #AP85751-1 50ug #AP85751-2 100ug #AP85751-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Melanocortin receptor 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	MC5-R
Accession No.	Uniprot ID:Q9MZV8
Uniprot	Q9MZV8
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MC5R

Note: This product is for in vitro research use only