

Mucin-1

Catalog No: #AP85755



Package Size: #AP85755-1 50ug #AP85755-2 100ug #AP85755-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mucin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	MUC-1,Episialin,CD227,Muc-1
Accession No.	Uniprot ID:Q02496
Uniprot	Q02496
GeneID	17829;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Muc1

Note: This product is for in vitro research use only