

Diphosphoinositol polyphosphate phosphohydrolase 3-beta

Catalog No: #AP85758

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP85758-1 50ug #AP85758-2 100ug #AP85758-3 1mg

Description

Product Name	Diphosphoinositol polyphosphate phosphohydrolase 3-beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	DIPP-3-beta,Diadenosine 5',5'''-P1,P6-hexaphosphate hydrolase 3-beta,Diadenosine hexaphosphate hydrolase (AMP-forming),Nucleoside diphosphate-linked moiety X motif 11,Nudix motif 11,DIPP3B
Accession No.	Uniprot ID:Q58CW0
Uniprot	Q58CW0
GeneID	616931;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NUDT11

Note: This product is for in vitro research use only