

## Metaxin-2

Catalog No: #AP85759



Package Size: #AP85759-1 50ug #AP85759-2 100ug #AP85759-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Metaxin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Mitochondrial outer membrane import complex protein 2,MNCb-0780
Accession No.	Uniprot ID:O88441
Uniprot	O88441
GeneID	53375;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Mtx2

Note: This product is for in vitro research use only