

Cell surface glycoprotein MUC18

Catalog No: #AP85760



Package Size: #AP85760-1 50ug #AP85760-2 100ug #AP85760-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cell surface glycoprotein MUC18
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Gicerin,Melanoma cell adhesion molecule,Melanoma-associated antigen MUC18,CD146,Muc18
Accession No.	Uniprot ID:Q9EPF2
Uniprot	Q9EPF2
GeneID	78967;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mcam

Note: This product is for in vitro research use only