

Mucin-like protein 1

Catalog No: #AP85761



Package Size: #AP85761-1 50ug #AP85761-2 100ug #AP85761-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mucin-like protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Protein BS106, Small breast epithelial mucin, SBEM, UNQ590/PRO1160
Accession No.	Uniprot ID:Q96DR8
Uniprot	Q96DR8
GeneID	118430;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: MUCL1

Note: This product is for in vitro research use only