

Transcriptional activator Myb

Catalog No: #AP85762



Package Size: #AP85762-1 50ug #AP85762-2 100ug #AP85762-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Transcriptional activator Myb
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Proto-oncogene c-Myb
Accession No.	Uniprot ID:P06876
Uniprot	P06876
GeneID	17863;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Myb

Note: This product is for in vitro research use only