

Mimecan

Catalog No: #AP85763



Package Size: #AP85763-1 50ug #AP85763-2 100ug #AP85763-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mimecan
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Osteoglycin, Osteoinductive factor, OIF, OIF, SLRR3A
Accession No.	Uniprot ID:P20774
Uniprot	P20774
GeneID	4969;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: OGN

Note: This product is for in vitro research use only