

Matrix metalloproteinase-28

Catalog No: #AP85764



Package Size: #AP85764-1 50ug #AP85764-2 100ug #AP85764-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Matrix metalloproteinase-28
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MMP-28,Epilysin,MMP25,UNQ1893/PRO4339
Accession No.	Uniprot ID:Q9H239
Uniprot	Q9H239
GeneID	79148;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MMP28

Note: This product is for in vitro research use only