

Matrix metalloproteinase-19

Catalog No: #AP85775



Package Size: #AP85775-1 50ug #AP85775-2 100ug #AP85775-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Matrix metalloproteinase-19
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MMP-19,Matrix metalloproteinase RASI,Matrix metalloproteinase-18,MMP-18,MMP18,RASI
Accession No.	Uniprot ID:Q99542
Uniprot	Q99542
GeneID	4327;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MMP19

Note: This product is for in vitro research use only