

Multidrug resistance-associated protein 6

Catalog No: #AP85776



Package Size: #AP85776-1 50ug #AP85776-2 100ug #AP85776-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Multidrug resistance-associated protein 6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ATP-binding cassette sub-family C member 6,Mrp6
Accession No.	Uniprot ID:Q9R1S7
Uniprot	Q9R1S7
GeneID	27421;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Abcc6

Note: This product is for in vitro research use only