

Serine/threonine-protein kinase mTOR

Catalog No: #AP85778



Package Size: #AP85778-1 50ug #AP85778-2 100ug #AP85778-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serine/threonine-protein kinase mTOR
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	FK506-binding protein 12-rapamycin complex-associated protein 1,FKBP12-rapamycin complex-associated protein,Mammalian target of rapamycin,mTOR,Mechanistic target of rapamycin,Rapamycin and FKBP12 target 1,Rapamycin target protein 1,FRAP,FRAP1,FRAP2,RAFT1,RAPT1
Accession No.	Uniprot ID:P42345
Uniprot	P42345
GeneID	2475;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MTOR

Note: This product is for in vitro research use only